

Lifting thread · Filler · Medical device · Beauty device · Skin care device

LINKING PEOPLE, SKILL, AND MACHINE IN BEAUTY INDUSTRY



Right to pursue Beauty

Derma Solution, as a global leader in the beauty industry such as lifting thread, fillers, beauty devices and others, is building a strong global infra (network) by exporting in 39 countries, and based on the various experiences and accumulated technologies, we strive to satisfy our customers.




MISSION

Ensuring everyone's "right to pursue Beauty"



VISION

A total solution company that connects people, technology, and machineries in the field of Beauty Industry.



VALUE

We develop our products, technologies, and machines, by considering the safety as our top priority for the "Healthy Beauty



Derma Solution

Derma Solution is a reliable global medical company chosen by famous Korean cosmetic surgeons, dermatologist, and celebrities

DERMA SOLUTION PRODUCT

WHITE ROSA LIFTING THREAD

ANOL'VE FILLER

MESO-INJECTION DEVICE

INTROFILL	16
INTROFILL EASY	17

SEMI-PERMANENT

MEDI TOUCH	19
MEDI TOUCH CLASSIC	19
MEDI TOUCH PREMIUM	20

MTS DEVICE

FN-1	22
FNS	22

BEAUTY DEVICE

BODY MASTER	24
POLKA RF	25
TRES	26
T-SONIC	27
COOL FIT	28
ENERGY UP	29
ZERO PEEL	30
MIUP HIFU	31
HIFU POINT	32
CELLUREON	33
MAGIC POWER	35
MAGIC CELL	36
MAGIC UP	37
MAGIC IMPACT · MAGIC POINT	38
BODY PLUS	38
DERMAWAVE	40
TRICELL · EXCEL LIPO	41
SKINPRO · SKINPROPLUS	42
MAGIC ICE · CRYOD	44
MAGIC PEEL · AQUA PLUS	45
ZIZIMI	46
SUNLIGHT MASTER ⁶³³	47



WHITE ROSA

lifting thread



THREAD MATERIALS



6Month~1Year

1Year~2Year

1Year~2Year

2Year~3Year




NEEDLE TYPE



ARROW

It's a two-way molding technique to fix strongly skin

WHITE ROSA ARROW is a two-way molding thread with molding technology. With its excellent tensile strength and strong skin fixation, the lifting effect and duration. It is characterized by strengthening and taking care of deeper and darker wrinkles. And after a certain amount of time, it is natural. It is made of PDO (melting thread) thread that decomposes, so the burden of side effects is low.



Two ways Molding room


Excellent Tensile strength

Strong Fixing my skin

LIFT KING

Power lifting effect of the special surface molding technique with sawtooth structure

WHITE ROSA LIFT KING is made by molding. It is a molding thread with a special surface structure with teeth. The sawtooth that got deeper in multiple directions is like a hook when pulling a string. It helps to hold skin and it's fixed. It has excellent fixation and strong lifting effects.



Saw-toothed structure Molding room

Strong Lifting effect


Stable Fixing power




KING COBRA

Effective lifting with a powerful 360° 6D spine

WHITE ROSA Lifting Thread, KING COBRA is a lifting thread containing bidirectional barbs in 360° 6D dimensionally. Stronger fixation and more elastic lifting effect can be expected as the KING COBRA is in the form of vertical barbs while existing COG threads were in the form of diagonal barbs.



 360° turn

Special Package


WHITE ROSA's special product

For all products except MONO, EYE, and SCREW,

18G needle + duoderm is included.



MONO




PDO

PCL

PLLA

GAUGE	NEEDLE	THREAD	USP
25G	25/50/60	30/70/90	4-0
26G	60/90	90/150	5-0
27G	38/50/60/90	50/70/90/150	5-0
29G	25/38/50/60	30/50/70/90	6-0
30G	12/25/38	12/30/50	7-0

SCREW




PDO

PCL

PLLA

GAUGE	NEEDLE	THREAD	USP
26G	60/90	60/150	5-0
27G	38/50/60/90	50/70/90/150	5-0
29G	25/38/50/60	30/50/70/90	6-0
30G	25/38	30/50	7-0

COG




PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
18G	100	160/170	2
19G	100	160/170	1
21G	60/90	120/160/210	2-0

KING COBRA



PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
18G	100	160/170	2/1
19G	100	160/170	1/0
21G	60/90	120/160	3-0

LIFT KING




PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
18G	100	160/170	3
19G	100	160/170	2/1
21G	60/90	120/160	0

ARROW




PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
18G	100	160/170	3
19G	100	160/170	2

BOMB




PDO

PCL

GAUGE	NEEDLE	THREAD	USP
23G	38/50/60	50/70/90	6-0

MESH



PDO

PCL

GAUGE	NEEDLE	THREAD	USP
19G	38/50/60	50/70/90	5-0

NOSE




PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
19G	38/50/60	50/70/90	2
21G	60	90	2-0

EYE




PDO

PCL

GAUGE	NEEDLE	THREAD	USP
30G	25	30	7-0

TWIN ARM



PDO

PCL

WPDO

GAUGE	NEEDLE	THREAD	USP
19G	100	400	2



ANOLVE

FILLER




01



ZERO
Irritation

02



Strong
Volume

03



Clinically
Tested

04



Long Lasting


05




Stable
shooting

ANOLVE meets all 5 conditions

INJECTION POINT




PRODUCT DETAIL



Product Classification	Single Use Medical Device	Single Use Medical Device	Single Use Medical Device
Packaging unit	1ml*1cc	1ml*1cc	1ml*1cc
HA Contents	24mg/ml	24mg/ml	24mg/ml
Lidocaine Content	0.3%	0.3%	0.3%
Needle	2*27g 1/2	2*27g 1/2	2*30g 1/2


ANOL'VE FILLER

Technology of ANOLVE Filler



Cross Linked Hyaluronic Acid
Natural combination through chain entanglement

Firm Gel Texture




With the firm gel texture which artificial crosslinking has been minimized, it provides naturalness and volumizing effect that has increased viscosity and cohesiveness.

Proven Safety



Through the 15 years of Clinical trials, ANOLVE FILLER is a brand which safety and efficacy are recognized, and its proven safety is the proud

Comparisons with others



Longevity of 4 times longer compared to others ANOLVE Filler is a standing out product in the global filler market and is loved all over world, including China, Vietnam, Japan, Canada, and others.

MESO-INJECTION DEVICE

Blooming your face



Moisture is basic! Shiny Glowing Skin!



INTROFILL

**Even improves skin tone, elasticity,
and fine wrinkles!**

Intro-fill is an automatic injection
Which injects directly into the dermis
evenly and accurately,
Using 9PIN needles of 31G/32G that
has less pain and maximized
the efficacy of the procedure
And its improvement of fine wrinkles,
whitening, moisturizing, elasticity, and
lifting effect is excellent



[Gun board (Handpiece) & needle]



[Advantages of Intro-fill]

▶ **Touchscreen**

Easy to use and easy to recognize UI
Auto injector function can be enabled

▶ **9 PIN multi-needle**

With 31G/32G 9 PIN multi-needles,
it's easy to mount or detach, and maximized
its effectiveness of the procedure

▶ **Accuracy**

By adjusting the length of the needle, it guarantees
the accurate depth of the procedure and it is light since
it's a gun type and wrist snap can be used as well

▶ **Less pain**

It uses 31G/32G 9PIN needle for injection designed with
3angles cutting to reduce pain. By using the 31G/32G 9 PIN
needles that was cut with 3 angles, it was designed for
easy injection and less pain

▶ **Various usages**

It can be used for intradermal hyaluronic acid injection,
fillers, toxins, obesity injections, and others.

▶ **Minimize Drug loss**

With accurate injector function and skin contact timing
calculation, the drug loss has been minimized.

Efficiency UP! Convenience UP! Intradermal HA Device

INTROFILL EASY

**With the convenient operation,
achieve Glowing Skin!**

Intro-fill easy is an automatic injecting system that injects hyaluronic acid and other various drugs by setting up different programs, which its 31G 5PIN needles maximize the procedure effect with less pain and its improvement of fine wrinkles, whitening, moisturizing, elasticity, and lifting effect is excellent.



[Intro-fill Easy product composition]



Main product body + Pump/Filter Other: Charger, 2 suction hose, additional pump module

[Intro-fill EASY advantages]

- ▶ **Minimize the pain in the procedure by installing the suction function**
- ▶ **As the Filter Module is replaceable, it is hygienic**
- ▶ **It used 31G 5PIN professional needle**
- ▶ **Any type of syringe can be attached**
- ▶ **It has lumirex fine-drug only program & all types of drugs can be injected**
- ▶ **Maximized convenience with wireless charger, high definition LCD, and ergonomic design.**

skin moisturizing

improves wrinkles

improves skin tone

increase elasticity

smooth skin

SEMI-PERMANENT

All about beauty



Excellent grip and simple design

MEDI TOUCH

More delicate procedure with digital handpiece that is slimmer, and more convenient to use.



1. adjustable Needle length (0.1-2.5)
2. Power Button / Speed control

- 01 Accurate, even and delicate touch is possible
- 02 Designed with advance mechanism to reduce pain
- 03 Compact size & ergonomic design
- 04 Performance UP with simple motion design (3-step Speed control)
- 05 Designed to minimize low noise/low vibration /operation temperature.
- 06 Uses Semi-permanent devices [indicator for medical usage]

Delicate Power control with 2 handpieces!

MEDI TOUCH CLASSIC

With slim, simple design and stable grip, delicate procedure can be done, and various procedures can be performed due to its needle processing technology and engineering design.



- 01 Compact size that is easy to move and carry
- 02 Induce absorption of medicines
- 03 Designed for low noise / low vibration / suitable operating temperature
- 04 9 modes / level 50~150 (5 increment) fine speed adjustable steps
- 05 Built-in RPM suitable for each procedure area, all control level
- 06 equipped with needle cartridge
- 07 2 handpieces can be connected

hypoallergenic indicator for medical use

MEDI TOUCH PREMIUM

Comfortable procedure with two handpieces!

Two handpieces can be used as indicator for medical usage and a drug absorption-inducing skin stimulator.



[Advantages of Meditouch Premium]

- 01** Two types of handpieces can be used (semi-permanent & MTS)
- 02** Easy and accurate needle length and RPM control
- 03** Easy to use
- 04** Digital built-in system with soft touch method
- 05** LCD Screen, for eye-comfort
- 06** Low noise/ low vibration/ low irritation
- 07** Simple and sophisticated design of Meditouch Premium

[Meditouch Premium Product Composition]

1. Main Body 2. Two types of handpieces (indicator for medical usage, drug absorption-inducing skin stimulator)
3. Handpiece connection cable 4. Adapter



※ Needles are separate product so its not included in the composition (separate purchase for needles)

MTS DEVICE

All about beauty



Fractional Needle therapy System

FN-1

A new paradigm in inducing drugs efficiently for the skin through microchannels

- 01** MTS devices that can be used in sophisticated procedures with the power from the strong motor
- 02** Micrometer method applied for needle exposure length (0.00mm~2.0mm precise adjustment)
- 03** Constant speed maintains in various procedure area
- 04** Little scratches only and small areas can be treated
- 05** 33G Pine-shaped needles is easy to permeate with less pain



hypoallergenic indicator for medical use

FNS

The microneedles (33G 9PIN needle) punctures and stimulates the skin to increase penetration of effective drugs into the skin through the puncture areas

- 01** MTS Handpiece that is easy to use and carry around
- 02** Adjustable speed step by step
- 03** Strong power and long lasting motor
- 04** Able to do detailed procedures in narrow areas
- 05** Able to do wide range of procedures in a short period of time
- 06** Comfortable grip and luxurious design



BEAUTY DEVICE

For beautiful skin



A collection of endermologie, high frequencies, suction, and cavitation.

BODY MASTER

Perfect on the outside, line management

high frequency, suction, and color therapy LED Bipolar, ultrasonic wave, fat layer only. Selection and management cavitation, continuous with the addition of a contraction and relaxation system, Energy transfer, complete line care



[BODY MASTER One-Stop System]

Care for details regardless of wide or narrow parts!

Take care of it with a variety of functions!

BODY MASTER is cellulite (orange peel), bumpy fold, forearm, It takes care of detailed areas such as thighs and sides. It helps to make wannabe body shape with various functions.


[BODY MASTER effect]



► Before

► After

[BODY MASTER, Systematic line management]



ROLL-RF-VACUUM

APP-S-RF-VACUUM

LED BIPOLAR

CAVITATION

Fractional Bipolar

POLKA RF

Care for both epidermis and dermis

Use dedicated tips where desired
Epidermal Regeneration Effects
Transferring energy to the dermis.
Implementing the regeneration
effect of high-frequency equipment




[POLKA RF Only Tools]

Heat stimulation to the skin layer
to shrink aging collagen.
Stimulates fibroblasts, helps induce
new collagen production.


[POLKA RF Principle]

Perfect care from epidermis to dermis.
Enough energy transferable!

Fine dots, high-frequency energy.
Divided and minimal damage
Maximize management effectiveness



[POLKA RF Advantages]



The epidermal regeneration effect of Fractional
and the high frequency dermis at once.
Helps induce collagen regeneration in the dermis.

**Safe management with
minimal management side effects**

CAVITATION + RF VACUUM

TRES

Total management of cavitation, high frequency, and thermal suction with one device
 Simultaneously manage cellulite decomposition and skin elasticity
 Manageable in a short period of time and selectively manage the fat layer



[CAVITATION + RF VACUUM]

Cellulite and Elasticity Management A combination for efficient total care!

Manage by selecting only the fat layer due to ultrasonic waves. A combination of cavitation and powerful suction high frequency. a sure before and after




CAVITATION



RF+VACUUM

[TRES Advantages]



Stimulate unnecessarily bloated fat by acting deep into the fat layer.
 Hard fat is softened and can be easily managed.

**Excellent effect on fat-rich
 abdomen, waist, hip, etc.**

A new paradigm for high-density ultrasound

T-SONIC

Elastic and moist skin with water droplets ultrasound


High-density ultrasound with wavelengths of 1, 3, 10, and 17MHz delivers energy to the entire skin layer.



[Various ultrasonic wavelengths unique to T-SONIC]

Four wavelengths deliver sufficient energy from shallow to deep layers of the skin.

With various wavelengths of 1, 3, 10 and 17MHz, T-SONIC can care for the deeper layers of the skin, such as the muscle layer and the fat layer, and even care for skin troubles in the shallow layer.




[Handpiece composition and utilization]



• **3, 10, 17MHz (Blue LED)**

Sensitive skin care
Various skin trouble care
Helps calm skin irritation



• **1, 3MHz (Red LED)**

Affects the muscle layer and the fat layer
to relax the muscles, lymphatic massage,
skin swelling removal, care before lifting

REDESIGN YOUR BODY

COOL FIT

COOL FIT's cooling system

For Cool Fit, after attaching the cool pad and handpiece to the area you want to be managed, it delivers cooling energy of minus 5-10 degrees to specific areas only.



Large Handpiece



Small Handpiece

[COOL FIT's cooling system]

**Stop the stressful diet anymore !
Start a comfortable diet with
COOL FIT's unique cooling system !**

COOL FIT delivers powerful cooling energy to the desired area such as the abdomen, sides, thighs, and arms to create the desired body shape.


[COOL FIT effect]



▶ Before

▶ After

[Changes in 90 days of using COOL FIT]



TO BALANCE YOUR BODY

ENERGY UP

Confirmed reliable care Rhythmic movement

In ENERGY UP, an energized manager receives a microcurrent using a living body current through a ground plate and delivers the microcurrent to a customer with a hand or handpiece.



[ENERGY UP Care Course]



The device transmits the microcurrent using the bio-current to the manager's feet through the ground plate.



The operator conducts a living body microcurrent to the ground plate.



The energized manager delivers the biological microcurrent to the customer by hand or handpiece.

[ENERGY UP effect]



▶ Before



▶ After

[Features of ENERGY UP]

Total care	Face and body total care possible by using biological microcurrent
Plates	Two managers can manage the upper and lower body with 4 ground plates
Handpiece	Two types of care systems that can be used with hand perfusion and pulse therapy systems.
Simple usage	Convenient management with simple usage

[PTS - 3 types of modes]

- MODE 1 : Care for local areas and lymph nodes
- MODE 2 : Systemic care for large areas, abdomen, etc.
- MODE 3 : Muscle care with multi-care mode

PEEL + E.P + P.T.S

ZERO PEEL

Finding the answer to my own skin troubles.

Separates dead skin cells with 25KHz-28KHz ultrasonic vibration, and removes dead skin cells in two types, dry and wet.

Contour·dead skin cells·skin troubles simultaneously care

- * Ultrasonic vibration of 25KHz-28KHz per second
- * With a hook-shaped tip, it is possible to care for curved areas of the face
- * Ampoule penetration deep into the skin
- * Facial line care with muscle stimulation



[Handpiece composition and utilization]



①	PEEL	Removes dead skin cells, blackheads, and whiteheads, penetrates ampoule into the skin
②	P.T.S E.P	Use for face narrow area such as face
③	P.T.S E.P	Use for body wide area such as body

[ZERO PEEL effect]



▶ Before

▶ After

[PTS effect]

Prevent muscle atrophy / relieve muscle spasm
 Facial lifting effect by stimulating the muscles of the face.
 Effective on reflex nerve stimulation such as body/foot floor

A new concept HIFU

MIUP HIFU

Reasons for flawless skin

MIUP HIFU is effective in skin elasticity and lifting care without damaging the skin surface.



[MIUP's high-intensity focused ultrasound]


Face and body management with high intensity focused ultrasound using the principle of a magnifying glass

Increased skin elasticity by applying high-intensity focused ultrasound to a specific area of the skin



*Ultrasound using the principle of a magnifying glass to treat only the desired area

[MIUP effect]



HIFU

POINT

Lifting without affecting everyday life

The HIFU system is a high intensity system that focused on ultrasound method of burning tissue using ultraviolet energy of 65-100°C, which occurs at the focus when ultrasound energy is collected in one place.


[Cartridge type]



[Product specification]

Type	D1.5	D3.0	D4.5
Frequency	7 MHz	4 MHz	4 MHz
Length	1.5mm	3.0mm	4.5mm
Life span	300,000 SHOT	300,000 SHOT	300,000 SHOT


[HIFU effect]



[Product standard]

Division	Standard
Rated Voltage	220V, 50/60Hz
Size	404.1 x 231.0 x 350.1(mm)
Weight	2.9kg

[Handpiece]




1. Operation indicator
2. Shot button
3. Cartridge

Heat Sensors

CELLUREON

New concept of decompression rolling system

Strong vibration generated by decompression capacity that is 15 times strong, makes rigid organ and muscles flexible. It is a new concept of decompression rolling system that facilitates blood circulation and promotes lymph circulation.



[CELLULION function]

Vacuum	a vacuum in the air inside the handpiece to reposition the fat and massage.
Suction	pulling caused by low pressure and vacuum. Technique and strength can be adjusted depends on degree of cup shape and suction.
Vibration	Generate heat inside the cell by repetitive movements periodically.
Acupressure	Make blood circulation and stimulate the nerve using nodules of the head.

[CELLULION body shaping]




It is a body management program that helps to manage cellulite by flexing rigid organ and muscles by improving circulation of blood and lymph by new concept of rolling system.

[CELLULION step 3]



[Thermotherapy]

Thermotherapy is a therapy that heals and improves the disease by maintaining some temperature of living organs or tissues in the range of 41°C - 45°C degrees within the limit to overcome the thermal effect while maintaining the physiological temperature control function.



► Before

► After



▶ 10.2 inch touchscreen



Flat type handpiece
- for use on back



Curved handpiece
- for use on arms and legs

Thermal
Vacuum



▶ Face handpiece




[Abdomen clinic once - 20min]



▶ Before



▶ After



▶ Before



▶ After

5 functions in one

MAGIC POWER

High frequency total care system

A system that combines the functions of RET, RF high frequency, Vacuum, LED, and Cryod into one.



[Handpiece name and function]



- 1 VACUUM+RF: For Body
- 2 VACUUM+RF: Small for body
- 3 VACUUM+RF: For Face
- 4 RET: For Body
- 5 RET: For Face
- 6 Scraping massage accessories
- 7 Scraping massage connector
- *Separately purchased
- 8 RET-only tips
- 9 Cryod Handpiece

[LED RF high frequency Vacuum / RET high frequency]



RET high frequency
 Helps generate bioenergy deep in the human body.
 Helps regenerate elastin collagen fibersz
 Helps in fat dissolution and body line management

MAGIC LINE+BIPOLAR

MAGIC CELL

Bipolar High Frequency + Vacuum + LED + Cryod

Efficient skin and body care with four functions of Bipolar High Frequency, Vacuum, LED, and Cryod in one device.



BIPOLAR effect

- 01 Cell function activation
- 02 Body fat and body size reduction
- 03 Cellulite emission
- 04 Pain relief
- 05 Relieve muscle atrophy


[Product series]



MAGIC SLIM

MAGIC LINE

[LED Bipolar suction]



Characteristics

Four types of suction mode available
1 to 10 RF level adjustment function
Independently available function

[MAGIC ALL]



- Vacuum
- RF high frequency
- LED

Total vacuum system

MAGIC UP

Bipolar High Frequency + Thermal + LED + Vacuum


Bipolar High Frequency + Thermal + LED + Vacuum
Efficient skin and body care with 4 functions in one device

MAGIC UP Advantages

- * Bipolar high frequency + LED + vacuum, thermal + vacuum bipolar high frequency, vacuum, etc. can be used in various modes.
- * No plate
- * Independently available function
- * Integrated effect of RF high frequency and vacuum
- * Total vacuum system



[MAGIC UP effect]



Powerful vibration and decompression function and recuperative energy heat work simultaneously to give physical stimulation to the target area for muscle relaxation, blood circulation and lymph circulation promotion, cellulite, lipolysis, elasticity, facial contouring and lifting effects.

[3 in MAGIC UP]

Vacuum	The powerful 650 mmHg of physical force stimulates blood and lymphatic flow; promoting metabolism and facilitating the evacuation of wastes.
LED	Body light Therapy = 630nm LED Face light Therapy = 430nm LED
RF	When a multipolar high-frequency current of 1.0MHz is energized in the human body, molecules of the tissue vibrate each other as the direction of current changes and cause deep heat. Improving skin elasticity, promoting fat dissolution and emission

[MAGIC UP accessories]



► MAGIC UP accessories



► Vacuum accessories



Decompression rolling system

MAGIC IMPACT

Care of stiff muscles through the action of strong and arious vibration modes and heat, and help to promote blood circulation and lymph circulation.

[MAGIC IMPACT effect]

- 01** Vacuum **02** Pain relief **03** Muscle relaxation
05 Improve immunity **05** Fatigue recovery



► Before

► After / once 25min care



Energy heat + Thermal vacuum

MAGIC POINT

Various vacuum accessories can care for various parts of the body.

Various applications are possible due to the thermal vacuum function.




► Vacuum accessories

CET RF

BODY PLUS

Care system using high-frequency energy

Designed to effectively deliver high-frequency energy to the skin deep into the epidermis and sebum fat layer. It helps to improve the treatment area quickly.




[CET RF product series]



BARON


DIARO

[CET RF high-frequency]



It is a method of contacting the anode with the insulation coated ceramic and the grounded conductive plate. There is active electron activity centering on the ceramics contacted with the skin, and deep heat is generated. Therefore, it is effective in managing the dermis and subcutaneous fat close to the skin tissue.

[CET RF high-frequency effect]



- Lipolysis · normalization of cell function
- Hormonal balance reulation
- Promote lymph circulation
- Promotes contraction and regeneration
- Stimulate collagen layer
- Promote blood circulation

BIPOLAR

DERMAWAVE

Third-generation of high-frequency system

The bipolar high-frequency system is a third-generation of high-frequency system that has been stripped from the existing high-frequency ceramics coating system.


It is an advanced system that enables strong high frequency and efficient in skin care along with body management without coating.



BIPOLAR effect

- 01 Cell function activation
- 02 Body fat and body size reduction
- 03 Cellulite emission
- 04 Pain relief
- 05 Relieve muscle atrophy

[Product features]



- 01 There is no spark that cause fragile ceramics and skin damage.
- 02 Easy to operate, fast and powerful.
- 03 It stimulates the energy of fat cells by deep heat.
- 04 We offer two handpieces, so it is easy to concentrate on each site.

[Clinical date]



▶ FACE

▶ FACE & BODY

▶ BODY

▶ BODY

CAVI+BIPOLAR


TRICELL

Cavitation Frequency 58KHz



▶ RF Handpiece

▶ Multi Handpiece




EXCEL LIPO

Find the answer to my own skin troubles

Ultrasound cavitation can be performed on all patients regardless of gender, and has excellent effects on obesity such as the fat-rich abdomen, waist, buttocks, and thighs.

[Egg burst simulation from cavitation power]



Comprehensive beauty machine

SKINPRO / SKINPROPLUS

All beauty care at once !

It is a multi-system with various skin care functions and various applicator for users can be provided for high treatment effect



- 01 Hygienic with single-use solution ampoule
- 02 Freely adjust solutions in aquafan
- 03 exfoliation, blackhead care, nutrition and moisturizing care at once
- 04 Reduced sebum, moist and clear skin tone
- 05 Simple and hygienic management

[SKINPRO]

Select only the desired applicator reconfigurable Minus Options System

**Request configuration of applicator without cryod - Possible*



- 1. Tongs pole · 2. Grounding pole · 3. Iontophoresis
- 4. Scrubber · 5. Ultrasound · 6. Spray

[SKINPROPLUS]



- 1. Aqua peel
- 2. Cryod
- 3. Grounding pole
- 4. Tongs pole
- 5. Iontophoresis
- 6. Scrubber
- 7. Ultrasound

[Velvet mask penetration · Spray]



▶ Velvet mask penetration tong



▶ Spray

Comprehensive beauty machine



- 3 time modes
- Sebum exfoliation
- 4 stage energy level adjustment



- 3 time modes
- Regeneration due to vitamin infiltration
- Anti whitening facial cleansing
- 4 stage energy level adjustment



- Crease regeneration by collagen penetration
- Wrinkle lifting
- 4 stage energy level adjustment



- Provides 2 programs
- Pulse function control
- 4 stage energy level adjustment



- Velvet mask penetration tong
- It can be applied safely and effectively by giving physical stimulation to the target area.
- 4 stage energy level adjustment



[SMASTER]

[For export SKINPROPLUS]

[Product series]



5-in-1 Skin Care

- 01 Skin scrubber mode
- 02 Iontophoresis mode
- 03 Lifting mode
- 04 Velvet mask penetration mode
- 05 Ultrasonic mode



BARON



DIARO



1. CRYO / 2. RF

Using the body's survival response

MAGIC ICE

CRYOD function and RF high frequency !

It is a principle to cleanse the body more cleanly by sending out unnecessary toxins and wastes in the body through the body's survival response.

MAGIC ICE system →

- 01 Powerful CRYOD therapy
- 02 Deep heat high frequency
- 03 Use of the two functions individually or together
- 04 Muscle relaxation
- 05 Cellulite reduction



- 1. Cold : -10°C / Warm : 55°C
- 2. Electroporation BODY
- 3. Electroporation FACE

CRYOD + HEATING + ELECTROPORATION

CRYOD

Electroporation

Electroporation applies an electrical pulse to the skin's surface to disrupt the skin's phospholipid membrane, temporarily creating a hole to penetrate the drugs needed to regenerate the skin.

Atopy relief Skin soothing Scalp care Pain relief



Supply plenty of moisture and
remove the waste of the face

MAGIC PEEL

Moisturize, remove face wastes!

AQUA PLUS PEEL uses a three step of special solution to inhale acne, blackhead, whitehead, exfoliation, folliculitis, dust, fine dust, etc. as the solution injected from the center of handpiece tip rapidly circulates based on the scientifically designed handpiece tip.




Electroporation



Skin care through water particles

AQUA PLUS

A system used for skin cleansing and cosmetic application by spraying moisture particles on the skin and scalp with powerful pressure



Far Infrared radiation steamer

ZIZIMI

If your body temperature goes up 1 degree, your immune system will increase 5 times.

Far infrared

It refers to infrared rays with a wavelength of 25mm or more. it is invisible because it has a longer wavelength than visible light. Also has a large heat effect and has strong penetration. Among these infrared rays, the longest wavelength is far-infrared.

[source : Doosan encyclopedia]



[Effect on heat action]



▶ Before




▶ After

it is effective to prevent various diseases such as aging, metabolism promotion, and chronic fatigue by activating the cell tissue by shaking the cell minutely 2,000 times a minute. When it comes into contact with moisture and protein molecules that make up the cell.

[Excellent Far-infrared emissivity]

[ZIZIMI Product standards]

Division	Standards	
Rated Voltage	220V, 60Hz	
Size	FLAT(mm)	510 × 80 × 360 (H)
	BEND(mm)	510 × 490 × 240 (H)
Weight	4.0 kg	



Excellent Far-infrared emissivity of 89.6% at 90°C

The Far-infrared emissivity rises in proportion to the temperature, therefore the emissivity at the actual operating temperature is important.

Whitening Management Program

SUNLIGHT MASTER⁶³³



Use German original 633nm wavelength lamp

- 01 Bluetooth function eliminates boredom during management
- 02 Promote collagen + elastin synthesis
- 03 Real 633nm wavelength lamp in Germany penetrates deep into the skin
- 04 Effect on dull skin, aging skin, acne, and troubled skin



► Before



► After



www.derma-solution.co.kr

T. +82. 31. 604. 1115

E. hello@derma-solution.co.kr

A. D-904, Hanam Techno Valley U1 Center, 947, Hanam-daero, Hanam-si, Gyeonggi-do, Republic of Korea