

EX
SQU MSS
MICRO SURFACE STIMULATION SYSTEM

PRECISION APPLICATION
PROFESSIONAL RESULTS

MSS is designed to support controlled, hygienic and consistent application of high-spec formulations.

DIRECT VIAL CONNECTION

- Reduced Product Loss
- Reduced Contamination Risk
- Immediate Application

SKIN APPLICATION GUIDE BY
NEEDLE LENGTH

0.25mm

PERSONAL / ENTRY-LEVEL SKIN CARE
Target layer: Epidermis

0.5 mm

AESTHETIC / ADVANCED SKIN CARE
Target layer: Sub- epidermal area

1.0 mm+

PROFESSIONAL / INTENSIVE SKIN PROGRAMS
Target layer: Dermal level

MSS TECHNICAL SPECIFICATIONS



CATEGORY	SPECIFICATION
PRODUCT NAME	MSS (MICRO SURFACE STIMULATION)
APPLICATION TYPE	NON-INVASIVE MICRO-SURFACE STIMULATION
CONNECTION TYPE	DIRECT VIAL CONNECTION
NEEDLE GAUGE	34G ULTRA-FINE
NEEDLE STRUCTURE	MICRO-NEEDLE ROLLER SYSTEM
INTENDED USE	PERSONAL / AESTHETIC
SOLUTION COMPATIBILITY	LOW-VISCOSITY, HIGH-SPEC FORMULATIONS
APPLICATION FOCUS	UNIFORM DISTRIBUTION & CONTROLLED STIMULATION



THE POWER OF
PLANT PDRN



COSMETICS
EXSOLU



NOVAMEDINTL



NOVAMEDINTL

NOVAMEDINTL
EXSOLU.NET
SALES.NOVAMED@GMAIL.COM

ROSE PDRN SERUM

DNA HYDRATION CARE

LIGHTWEIGHT. FRESH. BALANCED.

Rose-derived PDRN combined with Hyaluronic Acid helps maintain optimal moisture balance while supporting a fresh and comfortable skin condition throughout the day.

WHY ROSE PDRN?

Hydration is not simply about adding moisture. Healthy-looking skin begins with maintaining balance.

ROSE PDRN is designed to deliver lightweight hydration while helping skin feel calm, comfortable, and refreshed.

KEY COMPONENTS



ROSE PDRN

Plant-Derived DNA Conditioning Care



HYALURONIC ACID

Long-Lasting Moisture Retention



CENTELLA EV

Skin Comfort Support



VITAMIN B12

Healthy-Looking Skin Appearance

IDEAL FOR



HYDRATION




DAYTIME CARE



SKIN COMFORT



FRESH FINISH



GINSENG PDRN SERUM

DNA REGENERATION COMPLEX

DNA-BASED SKIN CONDITIONING

Formulated with Wild Ginseng-derived PDRN, Niacinamide, and Tranexamic Acid, this advanced serum helps support skin vitality, improve visible skin tone, and maintain a firmer, more refined-looking complexion while you sleep.

WHY GINSENG PDRN?

Healthy-looking skin is built through continuous recovery. Powered by Wild Ginseng-derived PDRN, this advanced formula helps support skin vitality, improve visible firmness, and promote a smoother, more refined-looking complexion overnight.

KEY COMPONENTS



GINSENG PDRN

Advanced DNA Skin Conditioning



NIACINAMIDE

Visible Tone Refining Care



TRANEXAMIC ACID

Complexion Clarity Support



HYALURONIC ACID

Moisture Retention Technology

IDEAL FOR



OVERNIGHT SKIN CARE



DULL & TIRED-LOOKING SKIN



FIRMNESS & ELASTICITY CARE



MATURE SKIN APPEARANCE



EXOSOME SKIN BOOSTER

ONE PLATFORM TOW TARGETS

ExoBene is an advanced exosome activation solution developed for both skin and scalp care, combining exosome-conditioned media, peptides, amino acids and conditioning factors.

EXOBENE-S SKIN-FOCUSED ACTIVATION

Supports:

- Skin hydration and moisture balance
- Skin barrier conditioning
- Surface vitality and texture improvement
- Professional skin activation routines

ExoBene-S combines exosome-conditioned media, peptide complexes, hyaluronic acid and NMF-supporting ingredients to help maintain healthy-looking skin.



EXOBENE-H SCALP-FOCUSED ACTIVATION

Supports:

- Scalp comfort and balance
- Follicle-supporting environment
- Scalp vitality
- Professional scalp conditioning routines

ExoBene-H combines exosome-conditioned media, biotin, inositol, amino acids, vitamins and peptides to support a healthy scalp environment.

